- DE LED Lichttherapiestift**
GB Light Therapy Pen
FR Stylo de luminothérapie avec technologie LED
IT Penna terapeutica luminosa a LED
ES Lápiz LED para fototerapia
PT Caneta luz de LED de terapia
NL Led lichttherapie pen
FI LED-valohoitokynä
SE LED ljusterapipenna
GR Μολύβι φωτοθεραπείας LED

Art. 85180

Gebrauchsanweisung
Instruction Manual
Mode d'emploi
Istruzioni per l'uso
Instrucciones de manejo

Bitte sorgfältig lesen!
Please read carefully!
A lire attentivement s.v.p.!
Da leggere con attenzione!
¡Por favor lea con cuidado!

Manual de instruções Por favor, ler cuidadosamente!
Gebruiksaanwijzing A.u.b. zorgvuldig lezen!
Käyttöohje Lue huolellisesti läpi!
Bruksanvisning Läsas noga!
Οδηγίες χρήσης Παρακαλούμε διαβάστε με προσοχή!

DE Gebrauchsanweisung	3
GB Instruction Manual	11
FR Mode d'emploi	19
IT Istruzioni per l'uso	27
ES Instrucciones de uso	35
PT Manual de instruções	43
NL Gebruiksaanwijzing	51
FI Käyttöohje	59
SE Bruksanvisning	67
GR Οδηγίες χρήσης	75

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- ❶ Blue Light Therapy
- ❷ On/Off and function switch
- ❸ Battery cover
- ❹ Red Light Therapy

IMPORTANT INFORMATION! RETAIN FOR FUTURE USE!

Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use. Should you give this device to another person, it is vital that you also pass on these instructions for use.

This instruction manual belongs to this device. It contains important information about starting up and operation. Read the instruction manual thoroughly. Non-observance of these instructions can result in serious injury or damage to the device.

WARNING

These warning notes must be observed to prevent any injury to the user.

CAUTION

These notes must be observed to prevent any damage to the device.

- Only use this unit according to the intended purpose and as directed in the instruction manual. The warranty is null and void if it is not used for the intended purpose.
- The device is not intended for commercial use.
- Check the device carefully for damage before each use. Do not use the device if you notice wear, damage or signs of improper use.
- Children must be supervised to ensure that they do not play with the device.
- Cleaning and user maintenance must not be performed by children unless supervised.

- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Do not attempt to repair the unit yourself in the event of any defects. Not only does this render the warranty null and void but you can also run the risk of serious hazards (fire, electric shock, injury). Only have repairs carried out by authorised service points.
- **Do not look directly into the light of the appliance to avoid damaging the eyes.**
- Consult your doctor first before using this unit if you have diabetes or other related illnesses.
- Pregnant women should note the important precautionary measures and the amount of stress they can take, consult your doctor first if necessary.
- Never use the device while pregnant or while breastfeeding.
- Do not treat any parts of the body that have swellings, burns, peeling skin, wounds or other sensitive points.
- The treatment should be pleasant. Stop at once and consult your doctor if you feel pain or the treatment seems unpleasant.
- Never use the device if you suffer from the following diseases or problems: Eczema, psoriasis, lesions or open wounds, infections, scars, ¹³

abrasions or skin burns, sunburn, redness or blistering, Diabetes or other metabolic diseases, skin cancer, potential malignoma, epilepsy, collagen disorders, herpes or blood clotting disorders, Immunodeficiency (e.g. HIV/Aids).

- The device may not be used after radiation or chemotherapy
- Use is also to be refrained from in case of photosensitivity or other light sensitivity (e.g. polymorphic light dermatosis, light urticaria, Lupus etc).

BATTERY SAFETY INSTRUCTIONS

- Do not disassemble batteries!
- Remove discharged batteries from the device immediately! Increased risk of leakage!
- Avoid contact with skin, eyes and mucous membranes! If battery acid comes in contact with any of this parts, rinse the affected area with copious amounts of fresh water and seek medical attention immediately!
- If a battery has been swallowed seek medical attention immediately!
- Insert the batteries correctly, observing the polarity!
- Always keep the battery compartment well closed!
- Remove the batteries from the device if it is not going to be used for an extended period!
- Keep batteries out of children's reach!
- Do not attempt to recharge these batteries! There is a danger of explosion!
- Do not short circuit! There is a danger of explosion!
- Do not throw into a fire! There is a danger of explosion!
- Keep unused batteries in their packaging away from metal objects in order to prevent short circuiting!

- 2.1 Scope of delivery and packaging** Please check first of all that the device is complete and is not damaged in any way. In case of doubt, do not use the appliance and contact your supplier or your service centre.
The following parts are included:
- 1 **medisana DC 300**
 - 1 Battery, 1.5 V, type AA
 - 1 instruction manual

The packaging can be reused or recycled. Please dispose properly of any packaging material no longer required. If you notice any transport damage during unpacking, please contact your supplier without delay.

WARNING

Please ensure that the polythene packing is kept away from the reach of children! Risk of suffocation!

- 2.2 Intended use** Targeted treatment of acne and impure skin with light therapy.
Blue light to improve acne. Red light, improves and alleviates the inflammation of the skin.

3 Operating

3.1 Functions

WARNING

Risk of injury! Do not use device longer than 5 minutes in one area in any mode and only use it twice a day.

Inserting/replacing batteries

Pull the battery cover **3** off the device.

Insert the battery (1.5 V, type: AA), please note the right polarity: Positive pole is facing towards the function button **2**. Reattach the battery cover.

Switch on the device

Switch the device on by pressing and holding the function button **2** a few seconds. On the thinner side of the device the blue therapy light is emitted **1**.

**continue the treatment**

Place the device close to your skin and let the blue light shine on the acne. Use the blue light twice a day for up to 5 minutes in one area.

If the device is switched on already, you may change the light colour by short press the button ②. The red therapy light is emitted on the thicker side of the device ④.

Place the device close to your skin and let the red light shine on inflammation of the skin, for example. Use the red light twice a day for up to 5 minutes in one area.

Switch off the device

The device is switched off automatically after 5 minutes. Press and hold function button ② for a few seconds to switch off manually.

Low battery warning

If battery is low, the device is switched off. Please change the battery.

- 4.1**
- Cleaning and maintenance**
- Switch off the device before cleaning.
 - Clean the unit with a soft cloth which has been lightly moistened with a mild soap. Never dip the unit into water for cleaning and ensure that no water gets into the unit at any time. Only use the unit again, when it is fully dry.
 - Store the device in a dry, cool place.

4.2 Disposal

This product must not be disposed together with the domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Consult your municipal authority or your dealer for information about disposal.

Battery disposal:

Do not dispose of old batteries with your household waste, but dispose of them at a battery collection station at a recycling site or in a shop.

4.3 Technical specifications

Name and model	: medisana light therapy pen DC 300
Power supply	: 1.5 V, 1 x 1.5 V, type: AA
Dimensions	: approx. 14.1 x 2.6 x 2.6 cm
Weight	: approx. 23 g
Article No.	: 85180
EAN number	: 40 15588 85180 3

In accordance with our policy of continual product improvement, we reserve the right for technical and design changes without notice.

The current version of this instruction manual can be found under
www.medisana.com

Warranty and repair terms

Please contact your supplier or the service centre in case of a claim under the warranty. If you have to return the unit, please enclose a copy of your receipt and state what the defect is. The following warranty terms apply: **1.** The warranty period for **medisana** products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice. **2.** Defects in material or workmanship will be corrected free of charge within the warranty period. **3.** Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts. **4.** The following is excluded under the warranty: **a.** All damage which has arisen due to improper treatment, e.g. non-observance of the user instructions. **b.** All damage which is due to repairs or tampering by the customer or unauthorised third parties. **c.** Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre. **d.** Accessories which are subject to normal wear and tear. **5.** Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

The service centre address is shown on the attached leaflet.

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medisana GmbH