DE / GB / FR / IT

MEDISANA®

DE Elektrischer Hornhautentferner CR 860
GB Electric callus remover CR 860

88570 05/2016

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GB Safety instructions

Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use.

Safety instructions

- The unit is intended for home use only.
Only use the unit for its intended purpose in accordance with the instructions for use.
This callus remover is intended to be used on feet, heels, hands and elbows only.

- Only replace with batteries of the same type, never use different types of batteries together or used batteries with new ones!
Insert the batteries correctly, observing the polarity!
Remove the batteries from the device if it is not going to be used for an extended period!

Device and controls

- Protective cap
Grinding attachment
Release button
ON/OFF button
Battery cover
Brush for cleaning

Items supplied

- 1 MEDISANA callus remover
2 grinding attachments
2 batteries (type AA)
1 brush
1 storage bag
1 instruction manual

The packaging can be reused or recycled. Please dispose properly of any packaging material no longer required.

Application

The MEDISANA callus remover is a new appliance which smoothes rough and calloused skin on the feet, hands and elbows.

Insert / replace batteries

Remove the battery compartment cover by lightly pressing the ribbed surface downwards and then pushing it down to remove it.

Starting up

The MEDISANA callus remover is supplied with a coarse and a fine grinding attachment.

Use

- Wash and clean the skin and dry thoroughly.
Make sure that a grinding attachment is inserted in the appliance.
Remove the protective cap, if the appliance is not ready to use.

Replaceable grinding attachments

The service life of a grinding attachment depends on the frequency and size of the area you are treating.

- Before you change a grinding attachment, make sure that the appliance is switched off.
Hold the grinding head with one hand and with the other one press the release button.

- Do not switch on the appliance if a grinding attachment has not been inserted.

Cleaning and care

- Clean the appliance after each use to ensure better performance.
Remove the protective cap from the grinding attachment.
Clean the callus remover using the brush.

Disposal

This product must not be disposed together with the domestic waste.

Specifications

Table with 2 columns: Name and model, Power supply, Speed, Operating conditions, Storage conditions, Dimensions, Weight, Article No, EAN Code, Accessories.

In accordance with our policy of continual product improvement, we reserve the right for technical and design changes without notice.

Warranty and repair terms

Please contact your supplier or the service centre in case of a claim under the warranty.

- The following warranty terms apply:
1. The warranty period for MEDISANA products is three years from date of purchase.
2. Defects in material or workmanship will be corrected free of charge within the warranty period.

MEDISANA AG Jaegersbergstraße 19 41468 NEUSS GERMANY

For service, accessories and spare parts, please contact: Murrays Health & Beauty School Lane Chandlers Ford Hampshire UK S0534YN

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Medisana AG